

Assessing the Impact of Integrating DermaSensor™ into Primary Care for Skin Cancer Detection and Management

Marcy Goisse, MD
Betty Thomas, MD

Montgomery Medical, Inc.
Offices in Smithfield and Brownsville, PA

Validated, FDA-cleared device to non-invasively test skin lesions for cancer^{1,2}

The DermaSensor Device

This non-invasive, handheld device captures and records spectral samples of skin lesions and, using a proprietary algorithm, classifies the given lesion's spectral properties against those of known malignant and benign lesions.

- Validated to detect melanoma, basal cell carcinoma, and squamous cell carcinoma^{1,2}
- Assesses skin lesions between 2.5 to 15 mm in diameter in <30 seconds^{1,2}
- Device classifies lesions as "Investigate Further" or "Monitor"
- An "Investigate Further" result is accompanied by a spectral score of 1–10, with higher scores having greater similarity to malignant lesions. A "Monitor" result is interpreted as low risk of malignancy^{1,2,3}
- Utility studies show the potential to reduce missed cancer assessments performed by PCPs by half³

Skin cancers diagnosed increased by 150% and dermatology referrals decreased by 13%

Implementation of DermaSensor into a Primary Care Practice

Drs. Goisse and Thomas reviewed skin lesion evaluation data and metrics for the 12 months prior to implementing DermaSensor into their practice compared to these same measures for the subsequent 3 months with use of the device. Data analyses were conducted using the practice's EHR data and billing records.

- Documented skin lesion evaluations increased from zero being the primary concern prior to introduction of the device to 34 primary skin lesion evaluations with use of the device
- Referrals to dermatology decreased 13%, from an average of 1.5 per month to 1.3 per month
- After device implementation, the clinic's rate of patients with eventually confirmed skin cancer diagnosis increased 150% among the reviewed population. The ultimate lesion diagnoses were provided via pathology following a biopsy procedure
- Given the 13% decrease in monthly dermatology referrals with the 150% increase in diagnosed skin cancers, the increase in skin cancer diagnoses were driven by an increase in PCP skin cancer sensitivity and an increase in their skin exam frequency. These findings suggest device use increases skin cancer detection while not increasing unnecessary referrals to dermatology



Documented skin lesion evaluations increased from zero in the prior 12 months to 34 in the 3-month period after introduction of the device

Insights into Billing Information in Visits With DermaSensor

The billing information analyzed found there were 34 skin lesion evaluations documented over the 3 months of device use with other billing insights as follows:

- These visits were billed as skin assessments alone (6) or in combination with other visit codes (28). Claims for these services were paid by various payers.
- Of the 34 skin check visits, 12 were reported using a modifier 25 with an annual wellness visit

The device use has been found to improve PCP confidence in studies and in this real-world evaluation³

Physicians Satisfaction and Confidence

Physicians' self-reported satisfaction with their skin exam capabilities increased from 40% and 50% prior to implementing DermaSensor into their practice to both being 90% with use of the device, using a scale of 0% being no satisfaction and 100% being complete satisfaction.

Their self-reported confidence in their skin exam capabilities similarly increased from 20% and 30% prior to device implementation to both being 100% with use of the device, using a scale of 0% being no confidence to 100% being complete confidence. A perceived anecdotal improvement in patient referral compliance was also reported.

Indications for Use: The DermaSensor™ device is indicated for use to evaluate skin lesions suggestive of melanoma, basal cell carcinoma, and/or squamous cell carcinoma in patients aged 40 and above to assist in the decision regarding referral of the patient to a dermatologist. The DermaSensor device should be used in conjunction with the totality of clinically relevant information from the clinical assessment, including visual analysis of the lesion, by physicians who are not already expertly trained in the diagnosis and management of skin cancer.

Disclosures: Funding for this white paper was provided by DermaSensor Inc. This white paper analyzed retrospective medical record and billing data. The information in this document does not constitute billing or coding advice, nor does it provide any guarantee of reimbursement.

References:

1. Hartman RI, Trepanowski N, Chang MS, et al. Multicenter prospective blinded melanoma detection study with a handheld elastic scattering spectroscopy device. J Am Acad Dermatol. 2024;15:24–31.
2. Merry SP, McCormick B, Nguyen VL, Chatha K, Croghan I, Leffell D. DERM-SUCCESS: Clinical Validation of an Elastic Scattering Spectroscopy (ESS) Device in Assisting Primary Care Physicians' Detection of Skin Cancer. J Clin Aesthet Dermatol 2023 Dec: 16(4 Suppl): s16
3. Seiverling EV, Agresta T, Cyr P, Caines L, Nguyen VL, Chatha K, Siegel DM. Clinical Utility of an Elastic Scattering Spectroscopy Device in Assisting Primary Care Physician's Detection of Skin Cancers. J Clin Aesthet Dermatol 2023 April: 16(4 Suppl): s16-17.